

MIRAI Medical
— *Quality of Life* —

product catalog.

www.miraii.co.id

daftar isi.

Agiography Sterile Set	
● Standard Pack - Paket 1.....	02
● Complete Pack - Paket 2.....	03
● Standard Pack + CAG Set - Paket 3.....	04
● Complete Pack + CAG Set - Paket 4.....	05
Angiographic Catheter.....	06
Guide Wire	
● Hydrophilic.....	07
● PTFE.....	07
Extension Tube.....	08
Pressure Bandage	
● Transradial I.....	09
● Transradial II.....	09
Manifold.....	10
Introducer Set	
● Transradial.....	11
● Femoral-Set.....	12
Inflation Device	
● Infaltion Device I.....	13
● Infaltion Device II.....	13
Y-Connector.....	13

angiography sterile set.

No #1 Accessories Cardiovascular in Indonesia

Memberikan perlindungan optimal dan kenyamanan bagi tim medis, dengan penerapan serta pengendalian terhadap infeksi secara efektif dan efisien. Angiography set juga dapat mencegah migrasi bakteri pasien infeksius melalui penggunaan drape dan gown dengan bahan yang berkualitas guna meningkatkan tingkat keamanan terhadap infeksi.



angioset paket 1.

Standard Pack



Angio Set.

- Surgical Gown Size L (SSMMS 45 GSM, PPE Laminasi 50 GSM)
- Angiography Drape 200 x 330 cm (SSMMS 45 GSM, PPPE Laminasi 50 GSM)
- Absorbent Towel 25 cm x 34 cm (Yichen)
- Body Cover 80 cm x 100 cm (Polyethylene 50 micron)
- Guardian Cover 80cm x 105 cm (Polyethylene 50 micron)
- Upper Image Cover 100 cm (Polyethylene 50 micron)
- Instrument Table Cover 140 cm x 180 cm (PPPE Laminasi 50 GSM)

angioset paket 2.

Complete Pack



- Angioset Standard Pack (1 pcs)

Disposable Set Miraii:

- Spuit Disposable Syringe with Needle Onemed 10 cc Luer Lock (2 pcs)
- Spuit Disposable Syringe without Needle Onemed 20 cc Luer Lock (1 pcs)
- Miraii Manifold (1 pcs)
- Miraii High Pressure 30 cm 1200 Psi (1 pcs)
- Miraii High Pressure 60 cm 1200 Psi (1 pcs)
- Nierbeken Bengkok Plastik 25 cm Biru Onemed (1 pcs)
- Bowl Round Plastic 150 ml (1 pcs)
- Bowl Round Plastic 500 ml (1 pcs)
- Handle Scaple Onemed No. 3 (1 pcs)
- Surgical Blade Sterile Serenity No. 11 (1 pcs)

angioset paket 3.

Standard Pack + CAG Set



- Angioset Standard Pack (1 pcs)

CAG Set Miraii:

- Miraii RA Band
- Miraii Femoral Introducer Sheath (6F/7F x 11 cm)
- Miraii Transradial Introducer Sheath (5F/6F x 7 cm)
- Miraii PTFE Guidewire J Tip (0.035"*150cm)
- Miraii Hydrophilic Coated Guidewire Angled Tip (0.038"*180cm)
- Miraii Angiography Catheter (6F, JL/JR, 3.5/4.0/ 5F TIG)

angioset paket 4.

Complete Pack + CAG Set




- Angioset Complete Pack (1 pcs)

CAG Set Miraii:

- Miraii RA Band
- Miraii Femoral Introducer Sheath (6F/7F x 11 cm)
- Miraii Transradial Introducer Sheath (5F/6F x 7 cm)
- Miraii PTFE Guidewire J Tip (0.035"*150cm)
- Miraii Hydrophilic Coated Guidewire Angled Tip (0.038"*180cm)
- Miraii Angiography Catheter (6F, JL/JR, 3.5/4.0/ 5F TIG)

angiographic. Catheter



Features:

- Soft tip reduces potential for vessel trauma
- Double braided stainless steel for 1:1 torque control
- Maximum pressure: 1200psi / 84 bar
- Maximum Guide wire: 0.038" / 0.97 mm

guide wire.

Hydrophilic




Features:

- High quality Niti core wire for better navigation and faster vessel selection
- TPU jacket with tungsten providing excellent visibility under X-ray
- Unique coating technology ensure a lubricant surface durable

guide wire.

PTFE



Features:

- High quality PTFE pre-coated stainless steel.
- Highly flexible and smoothly rounded wire tip.
- Special designed dispenser to achieve easy flush.
- J tip radius : 3.0 mm or 1.5mm.

extension tube.

Customizeable Length




Ref.No.	Length	Description	Ref.No.	Length	Description
PMT0101	30cm	500psi, No braided, F/M Luer Lock	PMT0101	30cm	1200psi, Braided, Rotating adapter
PMT0102	60cm		PMT0102	60cm	
PMT0103	90cm		PMT0103	90cm	
PMT0104	120cm		PMT0104	120cm	
PMT0105	150cm		PMT0105	150cm	
PMT0106	180cm		PMT0106	180cm	
PMT0801A	30cm	1200psi, No braided, F/M Luer Lock	PMT0901A	30cm	1200psi, No Braided, Rotating adapter
PMT0801B	50cm		PMT0901B	50cm	
PMT0801C	60cm		PMT0901C	60cm	
PMT0801D	90cm		PMT0901D	90cm	
PMT0801E	120cm		PMT0901E	120cm	
PMT0801F	150cm		PMT0901F	150cm	
PMT0801G	180cm		PMT0901G	180cm	
PMT0801H	200cm		PMT0901H	200cm	

pressure bandage.

Transradial I

- ◆ Medical grade silicone.
- ◆ Clear dial provides easy monitoring of puncture site.
- ◆ Silicone curved surface, less pressure achieving vascular hemostasia, more comfortable for patient.
- ◆ Pressure indication, convenient to adjust compression pressure for nursing staff.
- ◆ Comfortable band, available in three size: 18cm, 22cm, 26cm.



pressure bandage.

Transradial II



- ◆ Soft material fits well to patient's skin while easy to clean for the outer layer.
- ◆ Double balloon design to achieve better hemostatic effect.
- ◆ Puncture site visibility provided by transparent material.
- ◆ Comfortable band, available in three size: 24cm, 27cm.

manifold.

Pressure Rating Up to 500psi (34bar)



Ref.No.	Description
AM0201	2 ports ON with rotating adapter, right port orientation
AM0202	2 ports OFF with rotating adapter, right port orientation
AM0301	3 ports ON with rotating adapter, right port orientation
AM0302	3 ports OFF with rotating adapter, right port orientation
AM0303	3 ports ON with pressure line, right port orientation
AM0303A	3 ports OFF with pressure line, right port orientation
AM0307	3 ports ON with rotating adapter, left port orientation
AM0308	3 ports OFF with rotating adapter, left port orientation
AM0401	4 ports ON with rotating adapter, right port orientation
AM0402	4 ports OFF with rotating adapter, right port orientation
AM0209	Soft type 2 ports ON with rotating adapter, right port orientation
AM0210	Soft type 2 ports OFF with rotating adapter, right port orientation
AM0309	Soft type 3 ports ON with rotating adapter, right port orientation
AM0310	Soft type 3 ports OFF with rotating adapter, right port orientation

introducer set.

Transradial

- ◆ Longer taper ensure smooth transition between sheath and dilator.
- ◆ Atraumatic tip reduce the damage to the artery.
- ◆ Cross-cut haemostasis valve minimizing blood reflux and air aspiration.
- ◆ Stainless steel guidewire, hydrophilic guidewire with cannula needle available.
- ◆ Hydrophilic coated sheath available.




Ref.No.	Size	Length of Sheath	Guide wire	Needle	Hydrophilic Coating
TR01	5Fr	110cm	0.021"	21G/ 45mm	NO
TR02	6Fr	110cm			
TR02A	6Fr	160cm			
TR03	5Fr	110cm	0.021"	21G/ 45mm	NO
TR03A	5Fr	160cm			
TR04	6Fr	110cm			
TR04A	6Fr	160cm			
HTR03	5Fr	110cm	0.025"	20G/ Cannula needle	NO
HTR03A	5Fr	160cm			
HTR04	6Fr	110cm			
HTR04A	6Fr	160cm			
CTR03	5Fr	110cm	0.021"	21G/ 45mm	YES
CTR03A	5Fr	160cm			
CTR04	6Fr	110cm			
CTR04A	6Fr	160cm			
CHTR03	5Fr	110cm	0.025"	20G/ Cannula needle	YES
CHTR03A	5Fr	160cm			
CHTR04	6Fr	110cm			
CHTR04A	6Fr	160cm			

Note: The 0.025 inch guidewire is coated with Hydrophilic.

introducer set.

Set-Femoral

- Different color and number make easy identity for different French size.
- Atraumatic tip reduce the damage to the artery.
- Cross-cut haemostasis valve minimizing blood reflux and air aspiration.
- 0.035" or 0.038" Stainless steel guidewire, Ni-Ti guidewire available.
- Hydrophilic coated sheath available.



Ref.No.	Size	Length of Sheath	Guide wire	Needle	Hyrdophilic coating
IS05	4Fr	110cm	0.035"	18G/70mm	NO
IS01	5Fr	110cm	0.038"		
IS02	6Fr	110cm			
IS03	7Fr	110cm			
IS04	8Fr	110cm			
HIS05	4Fr	110cm	0.035"	18G/70mm	YES
HIS01	5Fr	110cm	0.038"		
HIS02	6Fr	110cm			
HIS03	7Fr	110cm			
HIS04	8Fr	161cm			

inflation.

Device



- Medium pressure(500psi) 3 ways stopcock with vented cap.
- 30cm (1200psi) extension tube.
- Volume: 20ml.

Ref.No.	Description
AI30	20ml, no leakage under 30ATM
AI40	20ml, no leakage under 40ATM

inflation.

Device II




- Medium pressure(500psi) 3 ways stopcock with vented cap.
- 30cm (1200psi) extension tube.
- Volume: 20ml.
- Ergonomic Button design, enable one-click, one-hand operation.
- Inflate 50% faster.

Ref.No.	Description
AI25	20ml, no leakage under 30ATM
AI25A	20ml, no leakage under 40ATM

y-connector.

Pack

- Largelumen: Max 9Fr(3.0mm) for various device compatibility.
- One-hand operation with 4 types: Rotating, Push-Click, Push-Pull and Push-Rotate.
- Clear polycarbonate body enable good observation of bleeding.
- Pressure injection allowance up to 400psi.



Ref.No.	Description
YCK111	Screw Y-Connector(1), Insertion tool(1), Torque Device(1)
YCK111A	Screw Y-Connector(1), Insertion tool(1), Torque Device(1), Extension tube with stopcocks(1)
YCK119B	Push-Click Y-Connector(1), Insertion tool(1), Torque Device(1)
YCK119A	Push-Click Y-Connector(1), Insertion tool(1), Torque Device(1), Extension tube with stopcocks(1)
YCK118B	Push-Pull Y-Connector(1), Insertion tool(1), Torque Device(1)
YCK118A	Push-Pull Y-Connector(1), Insertion tool(1), Torque Device(1), Extension tube with stopcocks(1)
YCK100C	Push-Rotate Y-Connector(1), Insertion tool(1), Torque Device(1)
YCK100F	Push-Rotate Y-Connector(1), Insertion tool(1), Torque Device(1), Extension tube with stopcocks(1)



quality of life.

For Further Information



Jakarta : Jl Munggang no. 86-87 Balekambang
Kramat Jati, Jakarta Timur, 13530



(021) 2204 7248
(+62) 821 8000 1102



info@miraii.co.id