



LumenCross® F2

OCT Imaging Catheter


Vivolight Official Partner **MEDITRANS** 

LumenCross® F2

Parameters

Model	Usable Length	Profile	Lens Position From Tip	Separate Lens Lumen
F2	1350mm	2.58Fr(0.86mm)	20mm	Yes

Excellent lesion crossing ability



Low-profile: 2.58 Fr




Hydrophobic coating




Enhanced out-layer of distal catheter

Outstanding capability to observe more distant artery



Short distance from Lens to Tip: 20mm

Mindray Official Partner **MEDITRANS** 

 **Jakarta** : Jl Munggang no. 86-87 Balekambang
Kramat Jati, Jakarta Timur, 13530

Surabaya : Rungkut Menanggal Harapan
Blok O No. 28, Surabaya 60293

Others : Bandung, Makassar, Medan, Balikpapan

 (+62) 21 2204 7248
(+62) 31 8786 2593
(+62) 821 8000 1102

 www.meditrans.id

 [meditrans.id](https://www.facebook.com/meditransid)

 info@meditrans.id

 [meditransid](https://www.instagram.com/meditransid)