


UltraFlyer

IVUS IMAGING CATHETER



Golden Eye in Coronary Vessel

- Broadband Frequency
- High Sensitivity
- Optimized Balance of Resolution & Penetration

MEDITRANS Contact Info



Jakarta : Jl Munggang no. 86-87 Balekambang
Kramat Jati, Jakarta Timur, 13530
Surabaya : Rungkut Menanggal Harapan
Blok O No. 28, Surabaya 60293
Others : Bandung, Makassar, Medan, Balikpapan



(+62) 21 2204 7248
(+62) 31 8786 2593
(+62) 821 8000 1102



www.meditrans.id



info@meditrans.id




[meditrans.id](https://www.facebook.com/meditrans.id)



[meditransid](https://www.instagram.com/meditransid)


≥ 6F Guide Catheter Compatible
0.014" Guide Wire Compatible

Official Partner **MEDITRANS** 



Golden Eye Transducer

- Central Frequency: 50 MHz
- Optimized Balance of Resolution & Penetration



High Resolution High Penetration

Smoothing Tips

- Exclusive smooth processing technology
- Lowering the risk of injury


Detachable Catheter

- More Solution for Sticking Problems



UltraFlyer

IVUS Catheter Introduction



Entry Profile	Compatible Guide Wire	Proximal Profile	Pullback Length	Rail Length	Working Length	Compatible Guide Catheter	Tip to Transducer Distance
2.7F 0.9mm	0.014"	3.3F 1.1mm	120mm	16mm	1370mm	6F	20mm