

UltraFlyer D01-65

IVUS CATHETER



Golden Eye in Coronary Vessel

- Broadband Frequency
- High Sensitivity
- Optimized Balance of Resolution & Penetration

Mindray Official Partner **MEDITRANS** 



Jakarta : Jl Munggang no. 86-87 Balekambang

Kramat Jati, Jakarta Timur, 13530

Surabaya : Rungkut Menanggal Harapan

Blok O No. 28, Surabaya 60293

Others : Bandung, Makassar, Medan, Balikpapan



(+62) 21 2204 7248

(+62) 31 8786 2593

(+62) 821 8000 1102



www.meditrans.id



info@meditrans.id




[meditrans.id](https://www.facebook.com/meditrans.id)



[meditransid](https://www.instagram.com/meditransid)


≥5F Guide Catheter Compatible
0.014" Guide Wire Compatible

Official Partner **MEDITRANS** 



Golden Eye Transducer

- Central Frequency: 65 MHz
- Optimized Balance of Resolution & Penetration



Smoothing Tips


- Exclusive smooth processing technology
- Lowering the risk of injury

Detachable Catheter

- More Solution for Sticking Problems

UltraFlyer[®]HD

IVUS Catheter Introduction



Tip Profile	Compatible Guide Wire	Proximal Profile	Pullback Length	Rail Length	Working Length	Compatible Guide Catheter	Tip to Transducer Distance
0.026"	0.014"	3.15F 1.05m	120mm	16mm	1370mm	5F	20mm