


Arrow® Thermodilution Catheters


When maneuverability and ease of insertion matter

Arrow Thermodilution Catheters from Teleflex help optimize patient care by measuring pressures and cardiac output, with an added capability to infuse fluids if needed. The thermodilution catheter enables easy insertion, handling, and management with the different firmness levels available for clinicians. Our torque-controlled catheters offer three body firmness levels – torque (firm), flexi torque (mid-range), and flotation (softest).

Four Lumen Thermodilution Catheter



Five Lumen Thermodilution Catheter



Ordering Information

Four Lumen

PRODUCT NO.	FRENCH SIZE	HEPARIN COATED	USEFUL LENGTH	NO. OF LUMENS	CATHETER BODY STYLE	RECOM'D ARROW INTRODUCER SIZE	BALLOON INFLATION CAPACITY	INFLATED BALLOON DIAMETER	MAXIMUM WIRE FOR INSERTION	CVP PROXIMAL LUMEN PORT DISTANCE FROM TIP	INFUSION LUMEN PORT DISTANCE FROM TIP	CE MARK
AI-07165	5 Fr.	No	80 cm	4	T	5 Fr.	1.0 cc	8 mm	0.018 in	15 cm	n/a	Yes
AI-07167	7 Fr.	No	110 cm	4	T	7 Fr.	1.5 cc	12 mm	0.025 in	29 cm	n/a	Yes
AI-07365-NHL	6 Fr.	No	110 cm	4	T	6 Fr.	1.0 cc	8 mm	0.021 in	29 cm	n/a	No
AI-07367	7 Fr.	No	110 cm	4	F	7 Fr.	1.5 cc	12 mm	0.025 in	29 cm	n/a	Yes
AI-07368-NHF	7 Fr.	No	110 cm	4	FL	7 Fr.	1.5 cc	12 mm	0.025 in	29 cm	n/a	No

Five Lumen with Infusion

PRODUCT NO.	FRENCH SIZE	HEPARIN COATED	USEFUL LENGTH	NO. OF LUMENS	CATHETER BODY STYLE	RECOM'D ARROW INTRODUCER SIZE	BALLOON INFLATION CAPACITY	INFLATED BALLOON DIAMETER	MAXIMUM WIRE FOR INSERTION	CVP PROXIMAL LUMEN PORT DISTANCE FROM TIP	INFUSION LUMEN PORT DISTANCE FROM TIP	CE MARK
AI-07377-NH	7.5 Fr.	No	110 cm	5	F	8 Fr.	1.5 cc	12 mm	0.025 in	29 cm	31 cm	No

For more information contact your local sales representative or customer service at +1.800.523.8446 or cs@teleflex.com

Note: Order quantity is one each.

Note: Catheter body style: T – Higher Torque Control, FL – Flexi Torque (mid-range), and F – Flotation (softest)

Each package includes: Thermistor catheter, control stroke syringe, non-vented dust cap, and extension line clamps.